



AKR1CL1 Polyclonal Antibody

Cat No: HR1AP6931

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | AKR1CL1 Polyclonal Antibody |
| Source | Rabbit |
| Applications | IF,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 37kDa |
| GenelD?Human? | AKR1CL1 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | AKR1CL1; Aldo-keto reductase family 1 member C-like protein 1 |
| Background | similarity:Belongs to the aldo/keto reductase family., |