

PEA3 Polyclonal Antibody

Cat No: HR1AP6293

For research use only

Overview

Product Name	PEA3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	ETV4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ETV4; E1AF; PEA3; ETS translocation variant 4; Adenovirus E1A enhancer-binding protein; E1A-F; Polyomavirus enhancer activator 3 homolog; Protein PEA3
Background	function:Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3'.,similarity:Belongs to the ETS family.,similarity:Contains 1 ETS DNA-binding domain.,