

Cleaved-Cathepsin D HC (L169) Polyclonal Antibody

Cat No: HR1AP2137

For research use only

Overview

Product Name	Cleaved-Cathepsin D HC (L169) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	CTSD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CTSD; CPSD; Cathepsin D
Background	cathepsin D(CTSD) Homo sapiens This gene encodes a member of the A1 family of peptidases. The encoded preproprotein is proteolytically processed to generate multiple protein products. These products include the cathepsin D light and heavy chains, which heterodimerize to form the mature enzyme. This enzyme exhibits pepsin-like activity and plays a role in protein turnover and in the proteolytic activation of hormones and growth factors. Mutations in this gene play a causal role in neuronal ceroid lipofuscinosis-10 and may be involved in the pathogenesis of several other diseases, including breast cancer and possibly Alzheimer's disease. [provided by RefSeq, Nov 2015],