

PAGE-3 Polyclonal Antibody

Cat No: HR1AP5931

For research use only

Overview

Product Name	PAGE-3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	PAGE3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PAGE3; GAGED1; Putative G antigen family D member 1; Prostate-associated gene 3 protein; PAGE-3
Background	PAGE family member 3(PAGE3) Homo sapiens This gene is a member of family of proteins that are expressed in a variety of tumors and in some fetal and reproductive tissues. Multiple alternatively spliced transcript variants have been observed. [provided by RefSeq, Jan 2015],