

p53AIP1 Polyclonal Antibody

Cat No: HR1AP8254

For research use only

Overview

Product Name	p53AIP1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	TP53AIP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TP53AIP1; p53-regulated apoptosis-inducing protein 1; p53AIP1
Background	tumor protein p53 regulated apoptosis inducing protein 1(TP53AIP1) Homo sapiens This gene is specifically expressed in the thymus, and encodes a protein that is localized to the mitochondrion. The expression of this gene is inducible by p53, and it is thought to play an important role in mediating p53-dependent apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Oct 2011],