

p53AIP1 Polyclonal Antibody Cat No: HR1AP8254

For research use only

Overview

| Product Name | p53AIP1 Polyclonal Antibody |
|--------------------------|--|
| Source | Rabbit |
| Applications | IHC-p,IF,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| lsotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | kDa |
| GeneID?Human? | TP53AIP1 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | TP53AIP1; p53-regulated apoptosis-inducing protein 1; p53AIP1 |
| Background | tumor protein p53 regulated apoptosis inducing protein 1(TP53AIP1) Homo sapiens This gene is specifically expressed in the thymus, and encodes a protein that is localized to the mitochondrion. The expression of this gene is inducible by p53, and it is thought to play an important role in mediating p53-dependent apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Oct 2011], |