

## KG19 Polyclonal Antibody Cat No: HR1AP5937

For research use only

## Overview

| Product Name          | KG19 Polyclonal Antibody  |
|-----------------------|---|
| Source                | Rabbit  |
| Applications          | WB,IHC-p,IF,ELISA   |
| Species Reactivity    | Human, Mouse  |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               | Dioeisa   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 27kDa   |
| GeneID?Human?         | BRI3BP  |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | BRI3BP; KG19; BRI3-binding protein; I3-binding protein; Cervical cancer 1 proto-oncogene-binding protein KG19; HCCRBP-1 |
| Background            | tissue specificity:Most abundantly expressed in brain, liver and kidney.,   |