

KG19 Polyclonal Antibody

Cat No: HR1AP5937

For research use only

Overview

Product Name	KG19 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	BRI3BP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BRI3BP; KG19; BRI3-binding protein; I3-binding protein; Cervical cancer 1 proto-oncogene-binding protein KG19; HCCRBP-1
Background	tissue specificity:Most abundantly expressed in brain, liver and kidney.,