

PEPC (Phospho Ser15) Polyclonal Antibody

Cat No: HR1AP9287

For research use only

Overview

Product Name	PEPC (Phospho Ser15) Polyclonal Antibody
Source	Rabbit
Applications	WB
Species Reactivity	Plants
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	105-110kDa
GeneID?Human?	PGC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PGC; Gastricsin; Pepsinogen C
Background	progastricsin(PGC) Homo sapiens This gene encodes an aspartic proteinase that belongs to the peptidase family A1. The encoded protein is a digestive enzyme that is produced in the stomach and constitutes a major component of the gastric mucosa. This protein is also secreted into the serum. This protein is synthesized as an inactive zymogen that includes a highly basic prosegment. This enzyme is converted into its active mature form at low pH by sequential cleavage of the prosegment that is carried out by the enzyme itself. Polymorphisms in this gene are associated with susceptibility to gastric cancers. Serum levels of this enzyme are used as a biomarker for certain gastric diseases including Helicobacter pylori related gastritis. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 1. [provided by RefSeq, Oct 2009],