

## c-Jun (7E12) Mouse mAb Cat No: HR1AM2449

For research use only

## Overview

Product Name	c-Jun (7E12) Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	43kDakDa
GeneID?Human?	3725
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Activator Protein 1 antibody, AP1 antibody, p39 antibody
Background	c-Jun is a member of the Jun Family containing c-Jun, JunB and JunD, and is a component of the transcription factor AP-1 (activator protein-1).