

## Gas1 Polyclonal Antibody Cat No: HR1AP6584

For research use only

## **Overview**

Product Name	Gas1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	GAS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GAS1; Growth arrest-specific protein 1; GAS-1
Background	growth arrest specific 1(GAS1) Homo sapiens Growth arrest-specific 1 plays a role in growth suppression. GAS1 blocks entry to S phase and prevents cycling of normal and transformed cells. Gas1 is a putative tumor suppressor gene. [provided by RefSeq, Jul 2008],