

## GAAP Polyclonal Antibody Cat No: HR1AP7662

For research use only

## Overview

Product Name	GAAP Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	TMBIM4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TMBIM4; GAAP; LFG4; CGI-119; Protein lifeguard 4; Golgi anti-apoptotic protein; Protein S1R; Transmembrane BAX inhibitor motif-containing protein 4; Z-protein
Background	similarity:Belongs to the BI1 family.,