

D55 Polyclonal Antibody

Cat No: HR1AP9037

For research use only

Overview

Product Name	D55 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	TPD52L3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TPD52L3; Tumor protein D55; hD55; Testis development protein NYD-SP25; Tumor protein D52-like 3
Background	tumor protein D52 like 3(TPD52L3) Homo sapiens This gene encodes a member of the tumor protein D52-like family of proteins. These proteins are characterized by an N-terminal coiled-coil motif that is used to form homo- and heteromeric complexes with other tumor protein D52-like proteins. The encoded protein may play a role in spermatogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009],