



CYB561D1 Polyclonal Antibody

Cat No: HR1AP6711

For research use only

Overview

Product Name	CYB561D1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	CYB561D1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CYB561D1; Cytochrome b561 domain-containing protein 1
Background	cofactor: Binds 2 heme groups non-covalently., similarity: Contains 1 cytochrome b561 domain.,