

DBC-1 Polyclonal Antibody Cat No: HR1AP8046

For research use only

Overview

Product Name	DBC-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	idioeisa
Clonality	
Concentration	1 mg/ml
Observed band	102kDa
GeneID?Human?	KIAA1967
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	KIAA1967; DBC1; DBIRD complex subunit KIAA1967; Deleted in breast cancer gene 1 protein; DBC-1; DBC.1; p30 DBC
Background	tissue specificity: Expressed ubiquitously in normal tissues. Expressed in 84 to 100% of neoplastic breast, lung, and colon tissues.,