

DAP-1 Polyclonal Antibody

Cat No: HR1AP6073

For research use only

Overview

Product Name	DAP-1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	DAP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DAP; DAP1; Death-associated protein 1; DAP-1
Background	death associated protein(DAP) Homo sapiens This gene encodes a basic, proline-rich, 15-kD protein. The protein acts as a positive mediator of programmed cell death that is induced by interferon-gamma. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, May 2014],