

# TH Polyclonal Antibody

**Cat No: HR1AP4712**

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | TH Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Applications          | WB,IHC-p,IF,ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Observed band         | 56kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| GeneID?Human?         | TH                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Alternative Names     | TH; TYH; Tyrosine 3-monooxygenase; Tyrosine 3-hydroxylase; TH                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Background            | tyrosine hydroxylase(TH) Homo sapiens The protein encoded by this gene is involved in the conversion of tyrosine to dopamine. It is the rate-limiting enzyme in the synthesis of catecholamines, hence plays a key role in the physiology of adrenergic neurons. Mutations in this gene have been associated with autosomal recessive Segawa syndrome. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Jul 2008], |