

# CSRNP3 Polyclonal Antibody

Cat No: HR1AP8720

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                      |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | CSRNP3 Polyclonal Antibody                                                                                                                                                                                                                                                           |
| Source                | Rabbit                                                                                                                                                                                                                                                                               |
| Applications          | WB,ELISA                                                                                                                                                                                                                                                                             |
| Species Reactivity    | Human,Mouse,Monkey                                                                                                                                                                                                                                                                   |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                      |
| Immunogen             |                                                                                                                                                                                                                                                                                      |
| Species               | Rabbit                                                                                                                                                                                                                                                                               |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                         |
| Isotype               |                                                                                                                                                                                                                                                                                      |
| Clonality             |                                                                                                                                                                                                                                                                                      |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                              |
| Observed band         | 64kDa                                                                                                                                                                                                                                                                                |
| GeneID?Human?         | CSRNP3                                                                                                                                                                                                                                                                               |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                      |
| Cellular localization |                                                                                                                                                                                                                                                                                      |
| Alternative Names     | CSRNP3; FAM130A2; TAIP2; Cysteine/serine-rich nuclear protein 3; CSRNP-3; Protein FAM130A2; TGF-beta-induced apoptosis protein 2; TAIP-2                                                                                                                                             |
| Background            | function: Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity. Plays a role in apoptosis.,sequence caution: Unusual initiator. The initiator methionine is coded by a non-canonical CTG leucine codon.,similarity: Belongs to the AXUD1 family., |