

Ezrin (phospho Tyr478) Polyclonal Antibody

Cat No: HR1AP8614

For research use only

Overview

Product Name	Ezrin (phospho Tyr478) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	EZR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81
Background	ezrin(EZR) Homo sapiens The cytoplasmic peripheral membrane protein encoded by this gene functions as a protein-tyrosine kinase substrate in microvilli. As a member of the ERM protein family, this protein serves as an intermediate between the plasma membrane and the actin cytoskeleton. This protein plays a key role in cell surface structure adhesion, migration and organization, and it has been implicated in various human cancers. A pseudogene located on chromosome 3 has been identified for this gene. Alternatively spliced variants have also been described for this gene. [provided by RefSeq, Jul 2008],