

Vimentin (Phospho-Tyr38) Polyclonal Antibody

Cat No: HR1AP9570

For research use only

Overview

Product Name	Vimentin (Phospho-Tyr38) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human:Y38,Mouse:Y38,Rat:Y38
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	VIM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	VIM
Background	vimentin(VIM) Homo sapiens This gene encodes a member of the intermediate filament family. Intermediate filamentents, along with microtubules and actin microfilaments, make up the cytoskeleton. The protein encoded by this gene is responsible for maintaining cell shape, integrity of the cytoplasm, and stabilizing cytoskeletal interactions. It is also involved in the immune response, and controls the transport of low-density lipoprotein (LDL)-derived cholesterol from a lysosome to the site of esterification. It functions as an organizer of a number of critical proteins involved in attachment, migration, and cell signaling. Mutations in this gene causes a dominant, pulverulent cataract.[provided by RefSeq, Jun 2009],