

DNA-PKCS (phospho Thr2638) Polyclonal Antibody Cat No: HR1AP7915

For research use only

Overview

Desident Name	
Product Name	DNA-PKCS (phospho Thr2638) Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	PRKDC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PRKDC; HYRC; HYRC1; DNA-dependent protein kinase catalytic subunit; DNA-PK catalytic subunit; DNA-PKcs; DNPK1; p460
Background	protein kinase, DNA-activated, catalytic polypeptide(PRKDC) Homo sapiens This gene encodes the catalytic subunit of the DNA- dependent protein kinase (DNA-PK). It functions with the Ku70/Ku80 heterodimer protein in DNA double strand break repair and recombination. The protein encoded is a member of the PI3/PI4-kinase family.[provided by RefSeq, Jul 2010],