




## Bcl-2 (phospho Thr74) Polyclonal Antibody

Cat No: HR1AP8132

For research use only

### Overview

|                       |  |
|-----------------------|--|
| Product Name          | Bcl-2 (phospho Thr74) Polyclonal Antibody  |
| Source                | Rabbit   |
| Applications          | IHC-p,ELISA  |
| Species Reactivity    | Human  |
| Recommended Dilutions |  |
| Immunogen             |  |
| Species               | Rabbit   |
| Storage               | -20°C/1 year   |
| Isotype               |  |
| Clonality             |  |
| Concentration         | 1 mg/ml  |
| Observed band         | kDa  |
| GenID?Human?          | BCL2   |
| Human Swiss-Prot No.  |  |
| Cellular localization |  |
| Alternative Names     | BCL2; Apoptosis regulator Bcl-2  |
| Background            | BCL2, apoptosis regulator(BCL2) Homo sapiens This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016], |



# bioelsa