

Caldesmon (phospho Ser789) Polyclonal Antibody

Cat No: HR1AP2630

For research use only

Overview

Product Name	Caldesmon (phospho Ser789) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	CALD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CALD1; CAD; CDM; Caldesmon; CDM
Background	caldesmon 1(CALD1) Homo sapiens This gene encodes a calmodulin- and actin-binding protein that plays an essential role in the regulation of smooth muscle and nonmuscle contraction. The conserved domain of this protein possesses the binding activities to Ca(2+)-calmodulin, actin, tropomyosin, myosin, and phospholipids. This protein is a potent inhibitor of the actin-tropomyosin activated myosin MgATPase, and serves as a mediating factor for Ca(2+)-dependent inhibition of smooth muscle contraction. Alternative splicing of this gene results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2008],