

Cleaved-Caspase-6 p11 (A194) Polyclonal Antibody

Cat No: HR1AP9555

For research use only

Overview

Product Name	Cleaved-Caspase-6 p11 (A194) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	CASP6 MCH2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	caspase 6, apoptosis-related cysteine peptidase
Background	caspase 6(CASP6) Homo sapiens This gene encodes a member of the cysteine-aspartic acid protease (caspase) family of enzymes. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic acid residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein is processed by caspases 7, 8 and 10, and is thought to function as a downstream enzyme in the caspase activation cascade. Alternative splicing of this gene results in multiple transcript variants that encode different isoforms. [provided by RefSeq, Oct 2015],