

NOXIN Polyclonal Antibody

Cat No: HR1AP5405

For research use only

Overview

Product Name	NOXIN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	85kDa
GeneID?Human?	NOXIN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NOXIN; C11orf82; Nitric oxide-inducible gene protein
Background	function:Acts as an anti-apoptotic factor and its absence increases cell death under normal and stress conditions. Can induce cell cycle arrest in the G1 or early S phase and this activity is independent of TP53/p53.,subcellular location:Accumulates in the nucleus in response to stress.,