

FAIM2 Polyclonal Antibody

Cat No: HR1AP4908

For research use only

Overview

Product Name	FAIM2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	FAIM2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FAIM2; KIAA0950; LFG; LFG2; NMP35; TMBIM2; Protein lifeguard 2; Fas apoptotic inhibitory molecule 2; Neural membrane protein 35; Transmembrane BAX inhibitor motif-containing protein 2
Background	function:Protects cells uniquely from Fas-induced apoptosis.,similarity:Belongs to the BI1 family.,subunit:Interacts with TNFRSF6.,tissue specificity:Widely expressed. Expressed at high levels in the brain especially in the hippocampus.,