




MUM1 Polyclonal Antibody

Cat No: HR1AP9582

For research use only

Overview

Product Name	MUM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	78kDa
GenID?Human?	MUM1 EXPAND1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	melanoma associated antigen (mutated) 1
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the MUM1 family.,



bioelsa