



TUSC2 Polyclonal Antibody

Cat No: HR1AP5748

For research use only

Overview

Product Name	TUSC2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GenelD?Human?	TUSC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TUSC2; C3orf11; FUS1; LGCC; PDAP2; Tumor suppressor candidate 2; Fusion 1 protein; Fus-1 protein; PDGFA-associated protein 2
Background	tumor suppressor candidate 2(TUSC2) Homo sapiens This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available. [provided by RefSeq, Jul 2008],