




TDE2L Polyclonal Antibody

Cat No: HR1AP4689

For research use only

Overview

Product Name	TDE2L Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GenID?Human?	SERINC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SERINC2; TDE2L; FKSG84; Serine incorporator 2; Tumor differentially expressed protein 2-like
Background	similarity:Belongs to the TDE1 family.,



bioelsa