

IRS-1 (phospho Ser1101) Polyclonal Antibody

Cat No: HR1AP2450

For research use only

Overview

Product Name	IRS-1 (phospho Ser1101) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	170kDa
GeneID?Human?	IRS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IRS1; Insulin receptor substrate 1; IRS-1
Background	insulin receptor substrate 1(IRS1) Homo sapiens This gene encodes a protein which is phosphorylated by insulin receptor tyrosine kinase. Mutations in this gene are associated with type II diabetes and susceptibility to insulin resistance. [provided by RefSeq, Nov 2009],