



PIG11 Polyclonal Antibody

Cat No: HR1AP4318

For research use only

Overview

Product Name	PIG11 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GenID?Human?	TP53I11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TP53I11; PIG11; Tumor protein p53-inducible protein 11; p53-induced gene 11 protein
Background	TP53I11 (Tumor Protein P53 Inducible Protein 11) is a Protein Coding gene.

