

PDRG1 Polyclonal Antibody

Cat No: HR1AP4282

For research use only

Overview

Product Name	PDRG1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	PDRG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PDRG1; C20orf126; PDRG; p53 and DNA damage-regulated protein 1
Background	function:May play a role in chaperone-mediated protein folding .,induction:By UV and repressed by TP53/p53.,similarity:Belongs to the prefoldin subunit beta family.,tissue specificity:Predominantly expressed in normal testis and exhibits reduced but detectable expression in other organs.,