

BCAS4 Polyclonal Antibody

Cat No: HR1AP7849

For research use only

Overview

Product Name	BCAS4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	BCAS4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BCAS4; Breast carcinoma-amplified sequence 4
Background	disease:A chromosomal aberration involving BCAS4 may be a cause of breast cancer. Translocation t(17;20)(q23;q13) with BCAS3.,similarity:Belongs to the cappuccino family.,tissue specificity:Brain, thymus, spleen, kidney and placenta. Overexpressed in most breast cancer cell lines.,