



SLC43A1 Polyclonal Antibody

Cat No: HR1AP4563

For research use only

Overview

Product Name	SLC43A1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	61kDa
GeneID?Human?	SLC43A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC43A1; LAT3; PB39; POV1; Large neutral amino acids transporter small subunit 3; L-type amino acid transporter 3; Prostate cancer overexpressed gene 1 protein; Solute carrier family 43 member 1
Background	solute carrier family 43 member 1(SLC43A1) Homo sapiens SLC43A1 belongs to the system L family of plasma membrane carrier proteins that transports large neutral amino acids (Babu et al., 2003 [PubMed 12930836]).[supplied by OMIM, Mar 2008],