

COP ?1 Polyclonal Antibody

Cat No: HR1AP3141

For research use only

Overview

Product Name	COP ?1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	COPZ1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	COPZ1; COPZ; CGI-120; HSPC181; Coatamer subunit zeta-1; Zeta-1-coat protein; Zeta-1 COP
Background	coatamer protein complex subunit zeta 1(COPZ1) Homo sapiens This gene encodes a subunit of the cytoplasmic coatamer protein complex, which is involved in autophagy and intracellular protein trafficking. The coatamer protein complex is comprised of seven subunits and functions as the coat protein of coat protein complex (COP)I-vesicles. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2012],