

# BET3 Polyclonal Antibody

Cat No: HR1AP2879

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                          |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | BET3 Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                 |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Applications          | WB,IHC-p,ELISA                                                                                                                                                                                                                                                                                                                                                                                                           |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                                                                          |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                             |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Observed band         | 22kDa                                                                                                                                                                                                                                                                                                                                                                                                                    |
| GeneID?Human?         | TRAPPC3                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Alternative Names     | TRAPPC3; BET3; CDABP0066; Trafficking protein particle complex subunit 3; BET3 homolog                                                                                                                                                                                                                                                                                                                                   |
| Background            | trafficking protein particle complex 3(TRAPPC3) Homo sapiens This gene encodes a component of the trafficking protein particle complex, which tethers transport vesicles to the cis-Golgi membrane. The encoded protein participates in the regulation of transport from the endoplasmic reticulum to the Golgi apparatus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012], |