



## ApoL2 Polyclonal Antibody

Cat No: HR1AP2794

For research use only

### Overview

|                       |                                                                                                                                                                                                                                                                                                                                                           |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | ApoL2 Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                 |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                    |
| Applications          | WB,IHC-p,IF,ELISA                                                                                                                                                                                                                                                                                                                                         |
| Species Reactivity    | Human                                                                                                                                                                                                                                                                                                                                                     |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                           |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                           |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                    |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                              |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                           |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                           |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                   |
| Observed band         | 37kDa                                                                                                                                                                                                                                                                                                                                                     |
| GeneID?Human?         | APOL2                                                                                                                                                                                                                                                                                                                                                     |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                           |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                           |
| Alternative Names     | APOL2; Apolipoprotein L2; Apolipoprotein L-II; ApoL-II                                                                                                                                                                                                                                                                                                    |
| Background            | apolipoprotein L2(APOL2) Homo sapiens This gene is a member of the apolipoprotein L gene family. The encoded protein is found in the cytoplasm, where it may affect the movement of lipids or allow the binding of lipids to organelles. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008], |