

NOL10 Polyclonal Antibody

Cat No: HR1AP4082

For research use only

Overview

Product Name	NOL10 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	NOL10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NOL10; Nucleolar protein 10
Background	similarity:Belongs to the WD repeat NOL10/ENP2 family.,similarity:Contains 5 WD repeats.,