

ABCA8 Polyclonal Antibody

Cat No: HR1AP2673

For research use only

Overview

Product Name	ABCA8 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	179kDa
GeneID?Human?	ABCA8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ABCA8; KIAA0822; ATP-binding cassette sub-family A member 8
Background	ATP binding cassette subfamily A member 8(ABCA8) Homo sapiens The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intracellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the ABC1 subfamily. Members of the ABC1 subfamily comprise the only major ABC subfamily found exclusively in multicellular eukaryotes. The encoded protein may regulate lipid metabolism and be involved in the formation and maintenance of myelin. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014],