



## TAUT Polyclonal Antibody

Cat No: HR1AP8330

For research use only

### Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | TAUT Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                           |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                             |
| Applications          | WB,IHC-p,ELISA                                                                                                                                                                                                                                                                                                                                                                                     |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                                                    |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                    |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                    |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                             |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                       |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                    |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                    |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                            |
| Observed band         | 70kDa                                                                                                                                                                                                                                                                                                                                                                                              |
| GeneID?Human?         | SLC6A6                                                                                                                                                                                                                                                                                                                                                                                             |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                    |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                    |
| Alternative Names     | SLC6A6; Sodium- and chloride-dependent taurine transporter; Solute carrier family 6 member 6                                                                                                                                                                                                                                                                                                       |
| Background            | solute carrier family 6 member 6(SLC6A6) Homo sapiens This gene encodes a multi-pass membrane protein that is a member of a family of sodium and chloride-ion dependent transporters. The encoded protein transports taurine and beta-alanine. There is a pseudogene for this gene on chromosome 21. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013], |