

ATP7B Polyclonal Antibody

Cat No: HR1AP7746

For research use only

Overview

Product Name	ATP7B Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	ATP7B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ATP7B; PWD; WC1; WND; Copper-transporting ATPase 2; Copper pump 2; Wilson disease-associated protein
Background	ATPase copper transporting beta(ATP7B) Homo sapiens This gene is a member of the P-type cation transport ATPase family and encodes a protein with several membrane-spanning domains, an ATPase consensus sequence, a hinge domain, a phosphorylation site, and at least 2 putative copper-binding sites. This protein functions as a monomer, exporting copper out of the cells, such as the efflux of hepatic copper into the bile. Alternate transcriptional splice variants, encoding different isoforms with distinct cellular localizations, have been characterized. Mutations in this gene have been associated with Wilson disease (WD). [provided by RefSeq, Jul 2008],