

MDM1 Polyclonal Antibody

Cat No: HR1AP7970

For research use only

Overview

Product Name	MDM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	MDM1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MDM1; Nuclear protein MDM1
Background	Mdm1 nuclear protein(MDM1) Homo sapiens This gene encodes a nuclear protein similar to the mouse double minute 1 protein. The mouse gene is located in double minute (DM) chromatin particles, is amplified in the mouse transformed 3T3 cell line, and the encoded protein is able to bind to p53. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2011],