



Hrs (phospho Tyr334) Polyclonal Antibody

Cat No: HR1AP9079

For research use only

Overview

Product Name	Hrs (phospho Tyr334) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	86kDa
GeneID?Human?	HGS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HGS; HRS; Hepatocyte growth factor-regulated tyrosine kinase substrate; Hrs; Protein pp110
Background	hepatocyte growth factor-regulated tyrosine kinase substrate(HGS) Homo sapiens The protein encoded by this gene regulates endosomal sorting and plays a critical role in the recycling and degradation of membrane receptors. The encoded protein sorts monoubiquitinated membrane proteins into the multivesicular body, targeting these proteins for lysosome-dependent degradation. [provided by RefSeq, Dec 2010].