

ABCA6 Polyclonal Antibody

Cat No: HR1AP2672

For research use only

Overview

Product Name	ABCA6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	185kDa
GeneID?Human?	ABCA6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ABCA6; ATP-binding cassette sub-family A member 6
Background	ATP binding cassette subfamily A member 6(ABCA6) Homo sapiens The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intracellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, and White). This encoded protein is a member of the ABC1 subfamily. Members of the ABC1 subfamily comprise the only major ABC subfamily found exclusively in multicellular eukaryotes. This gene is clustered among 4 other ABC1 family members on 17q24 and may play a role in macrophage lipid homeostasis. [provided by RefSeq, Jul 2008],