



DBP Polyclonal Antibody

Cat No: HR1AP3260

For research use only

Overview

Product Name	DBP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	GC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GC; Vitamin D-binding protein; DBP; VDB; Gc-globulin; Group-specific component
Background	GC, vitamin D binding protein(GC) Homo sapiens The protein encoded by this gene belongs to the albumin gene family. It is a multifunctional protein found in plasma, ascitic fluid, cerebrospinal fluid and on the surface of many cell types. It binds to vitamin D and its plasma metabolites and transports them to target tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Feb 2011],