

SLC6A14 Polyclonal Antibody

Cat No: HR1AP4943

For research use only

Overview

Product Name	SLC6A14 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	72kDa
GeneID?Human?	SLC6A14
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC6A14; Sodium- and chloride-dependent neutral and basic amino acid transporter B ⁰ +; Amino acid transporter ATB ⁰ +; Solute carrier family 6 member 14
Background	solute carrier family 6 member 14(SLC6A14) Homo sapiens This gene encodes a member of the solute carrier family 6. Members of this family are sodium and chloride dependent neurotransmitter transporters. The encoded protein transports both neutral and cationic amino acids. This protein may also function as a beta-alanine carrier. Mutations in this gene may be associated with X-linked obesity. A pseudogene of this gene is found on chromosome X.[provided by RefSeq, May 2010],