

## **HORMAD1 Polyclonal Antibody**

Cat No: HR1AP3666

For research use only

## **Overview**

Product Name	HORMAD1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
sotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	HORMAD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HORMAD1; NOHMA; HORMA domain-containing protein 1; Cancer/testis antigen 46; CT46; Newborn ovary HORMA protein
Background	HORMA domain containing 1(HORMAD1) Homo sapiens This gene encodes a HORMA domain-containing protein. HORMA domains are involved in chromatin binding and play a role in cell cycle regulation. The encoded protein may play a role in meiosis, and expression of this gene is a potential marker for cancer. A pseudogene of this gene is located on the long arm of chromosome 6. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2010],