

LDOC1L Polyclonal Antibody

Cat No: HR1AP3819

For research use only

Overview

Product Name	LDOC1L Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human, Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	LDOC1L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LDOC1L; MAR6; MART6; Protein LDOC1L; Leucine zipper protein down-regulated in cancer cells-like; Mammalian retrotransposon-derived protein 6
Background	similarity:Belongs to the LDOC1 family.,