

# Galectin-4 Polyclonal Antibody

Cat No: HR1AP3505

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                  |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Galectin-4 Polyclonal Antibody                                                                                                                                                                                                                                                                                                   |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                           |
| Applications          | WB,IHC-p,ELISA                                                                                                                                                                                                                                                                                                                   |
| Species Reactivity    | Human                                                                                                                                                                                                                                                                                                                            |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                  |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                  |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                           |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                     |
| Isotype               |                                                                                                                                                                                                                                                                                                                                  |
| Clonality             |                                                                                                                                                                                                                                                                                                                                  |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                          |
| Observed band         | 36kDa                                                                                                                                                                                                                                                                                                                            |
| GeneID?Human?         | LGALS4                                                                                                                                                                                                                                                                                                                           |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                  |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                  |
| Alternative Names     | LGALS4; Galectin-4; Gal-4; Antigen NY-CO-27; L-36 lactose-binding protein; L36LBP; Lactose-binding lectin 4                                                                                                                                                                                                                      |
| Background            | galectin 4(LGALS4) Homo sapiens The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. The expression of this gene is restricted to small intestine, colon, and rectum, and it is underexpressed in colorectal cancer. [provided by RefSeq, Jul 2008]. |