

FASTKD2 Polyclonal Antibody

Cat No: HR1AP3442

For research use only

Overview

Product Name	FASTKD2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	81kDa
GeneID?Human?	FASTKD2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FASTKD2; KIAA0971; FAST kinase domain-containing protein 2
Background	FAST kinase domains 2(FASTKD2) Homo sapiens This gene encodes a protein that is localized in the mitochondrial inner compartment and that may play a role in mitochondrial apoptosis. Nonsense mutations have been reported to result in cytochrome c oxidase deficiency. [provided by RefSeq, Oct 2008],